
Instructions for use Mofixx Trocar

EN/NL/DE/
HU/SK/SL

Manufacturer:

Mofixx B.V.
Josink Esweg 1
7545 PN Enschede, the Netherlands

P.O.Box 265
7500 AG Enschede, the Netherlands

SRN: NL-MF-000008385

Telephone: +31(0)53 480 39 27

Website: www.mofixx.com

Email: info@mofixx.com

Service department of your local sales agent:

See website: www.mofixx.com

CE 0344 ©Copyright 2023-12-06; Version: 9.00

Mofixx B.V. owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mofixx is strictly forbidden.

Instructions for Use Mofix Trocar

EN

1 DEVICE DESCRIPTION

The Mofix trocar (see figure on front page) consists of four parts:

1=obturator; 2=trocar; 3=reusable metal ball clamp; 4= clamp base.

The Mofix Trocars is designed to provide a port of entry for laparoscopic cameras and instruments which are used during laparoscopic (keyhole) surgery, which is a procedure designed to cause minimal invasion to a patient requiring abdominal surgery.

The Mofix trocar is designed to be used in combination with the Mofix System. The design of the Mofix trocar deviates from generally available trocars, by the addition of a sphere that serves as the fixation point for the ball clamp of the Mofix System and enables a laparoscopic camera or instruments to be fixed in a desired position.

The clear, ribbed cannula provides for instrument visualisation and abdominal wall retention. The universal seal accommodates 10-10.3 mm instruments whilst maintaining pneumoperitoneum. A stopcock valve is compatible with standard luer lock fittings and provides attachment for gas insufflation and desufflation.

Basic UDI-DI: 87202997147Z1101- TROCAR9E

The trocar comes in four variants:

Catalogue number	UDI-DI	Shaft length	For camera or instrument diameter	Sphere inner diameter	Sphere colour
Z1101-1300	08720299714717	64.4 mm (short)	10 mm	10.2 ± 0.1 mm	Purple
Z1101-1400	08720299714724	96 mm (long)	10 mm	10.2 ± 0.1 mm	Purple
Z1101-1310S	08720299714731	64.4 mm (short)	10.3 mm	10.5 ± 0.1 mm	Colourless
Z1101-1410S	08720299714748	96 mm (long)	10.3 mm	10.5 ± 0.1 mm	Colourless

The Mofix Trocar is classified as a medical device according to the European Directive for Medical Devices 93/42/EEC.

2 INTENDED USE

The Mofix trocar provides a port of access to the abdomen for an endoscope or laparoscopic surgical instrument.

The Mofix trocar can be used with or without the Mofix holding system.

The Mofix trocar is intended for single use.

3 INDICATIONS FOR USE

Patients undergoing abdominal, urological and/or gynaecological minimally invasive surgery procedures.

4 CONTRAINDICATIONS

The Mofix disposable trocar is not intended for use when endoscopic techniques are contraindicated.

Instructions for Use Mofixx Trocar

EN

There is no contraindication that has direct connection with the Mofixx disposable trocar till now. Physicians should determine whether this product can be used for the patient's particular situation.

5 ACCESSORIES AND OTHER DEVICES TO BE USED WITH THE DEVICE

Following single use accessories are packed with the Mofixx Trocar

- Obturator with a bladeless configuration
- Clamp base, only required when used in combination with the Mofixx System

Following can be used with the Mofixx Trocar

- Mofixx System; a laparoscopic camera and instrument holder
- Laparoscopic camera and instrument

6 WARNINGS

- Inappropriate use of this device may result in injury to blood vessels and to the internal organs of the patient or may cause injury to the physician. Use extreme caution when inserting this device into the patient.
- The Mofixx Trocar, including obturator are intended for single patient use ONLY. Never reuse a single-use device. Re-use can cause cross-infection between patients.
- Do not attempt to resterilise.
- The Mofixx Trocar is supplied nonpyrogenic and sterile using an Ethylene Oxide (EO) sterilisation method. Never use the Mofixx Trocar after the expiry date indicated on the label, unintentionally opened before use, or with a damaged sterile barrier packaging.
- Careful handling of this device is necessary to avoid damage or breakage. Always inspect the device prior to use to ensure device integrity and do not use the device if any damage is identified. Return damaged devices along with its packaging to Mofixx for evaluation.
- Do not attempt insertion of the Mofixx trocar until the primary site and recommended pneumoperitoneum are established.
- Direct the trocar away from major vessels and other anatomic structures. Puncture of internal organs or major blood vessels may result in life-threatening complications.
- Do not use excessive force or in a manner not consistent with normal instrumentation use. This may compromise its functionality and safety.
- The device is intended for use with the Mofixx holding system. Never use the Mofixx Trocar in combination with clamping devices other than the Mofixx ball clamp and System.
- Keep dry and away from the sun.

7 PRECAUTIONS

- Use of the Mofixx Trocar is restricted to surgeons who have had relevant and adequate training on the use of trocars, are experienced in performing abdominal, urological, and gynaecological minimally invasive surgery procedures, and are familiar with the conceivable complications that may occur at any time during or after the procedure. If further information is required, consult medical literature regarding use, techniques, complications, and hazards prior to the use of this device.
- Prior to use, read and follow the instructions of all the products to be used during the procedure.
- Before opening the sterile barrier, verify that labelled dimensions are suitable for the intended laparoscopic camera or instrument.
- Prior to use, ensure compatibility with all devices intended to be used with the Mofixx trocar during the surgical procedure. Check if the laparoscopic camera or instrument intended to be used moves smoothly through the trocar when the release button of the Mofixx switch is pressed, or is properly fixated when the release button of the Mofixx switch is not pressed.
- Select and follow a suitable endoscopic and/or laparoscopic protocol.
- The trocar is delivered with the port stopcock in open position. Close the stopcock prior insertion to prevent de-sufflation.

mofixx

Instructions for Use Mofixx Trocar

EN

- To prevent the premature dislodging of the trocar, use care when removing instruments.
- Following use in a surgical procedure the product may pose a biological hazard. Handle and dispose the product according to local applied regulations for bio-hazardous waste.

8 DIRECTIONS FOR USE

- The trocar and obturator and clamp base are packaged unassembled. Assemble the trocar by inserting the obturator into the trocar cannula until they match together.
- Make a skin incision appropriate for the trocar diameter at the site where the trocar will be introduced. A greater insertion force will be required if the incision is too small, and if the incision is too large, the port may become unstable.
- Ensure there is no blood vessel or grume existing in the incision area. Insert the obturator and trocar according to the texture level of tissue through the incision by applying continuous downward pressure until the abdominal wall has been completely penetrated.
- Keep organs out of reach of trocar penetration by ensuring proper positioning of the patient's body.
- Direct the trocar away from major vessels and other anatomic structures.
- To insufflate, attach a gas line to the trocar port and open its valve.
- For the second and additional punctures of the trocar into the abdominal or thoracic cavity, inspect the tip of the trocar visually by monitor.
- Remove the obturator and insert the appropriately sized camera or instrument.
- Minimize device exchange.
- Leave the laparoscope in place during de-sufflation and removal of the trocar. Exteriorization of the cavity contents can occur if the laparoscope is first pulled from the trocar.
- Exhaust maximum gas before removing the trocar.
- Remove the trocar as efficiently as possible by using continuous upward force while gentle rotating the trocar in alternating clockwise and counter clockwise directions until the cannula is completely visible.
- After removing the trocar from the cavity, take appropriate steps to achieve haemostasis as needed.

9 COMPLICATIONS

Possible complications include, but are not limited to the following:

- Vascular injuries
- Visceral injuries
- Post-operative trocar site hernia (PTSH)
- Surgical site infections
- Allergic or toxic reaction
- Aseptic inflammatory reactions of fat layer
- Extraperitoneal /subcutaneous emphysema
- Extraperitoneal insufflation
- Haemorrhage
- Incision pain
- Failed entry

Any serious incident (incident that directly or indirectly led, might have led or might lead to death, temporary or permanent health deterioration of a patient, user or any other person or a serious public health risk) that has occurred related to the Mofixx system, should be reported to the manufacturer and / or distributor and to the competent authority in the Member State in which the user is located.

10 EXPLANATION OF SYMBOLS

One or more of the following markings may appear on the device

Symbol	Explanation
	CE mark including identification notified body
	Medical device
	Manufacturer
	Catalogue number
	Batch code
	Unique Device Identification
	Use by date
	Do not use if package is damaged and consult instructions for use
	Caution
	Consult instructions for use or consult electronic instructions for use
	Do not re-use
	Do not re-sterilise
	Sterilized using ethylene oxide
	Single sterile barrier system
	Packaging unit

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