



mofixx

---

Execute minimal invasive surgery with improved quality and improved cost efficiency.

- ✓ Better ergonomics in OR
- ✓ cost reduction per procedure
- ✓ more stable images
- ✓ compact
- ✓ "solo surgery" possible
- ✓ versatile (for endoscope & instrument)
- ✓ "drag and drop"



## Components

### Single use:

- ✓ Trocar
- ✓ Cover
- ✓ Switch

### Re-useable:

- ✓ Adapter (with actuator)
- ✓ Arm
- ✓ Power cable
- ✓ Electronics box
- ✓ Medical power supply
- ✓ Clamp

## Advantages:

Mofix is a product that enables surgeons to execute MIS procedures with improved quality and improved cost-efficiency. The benefits of using Mofix are the improved working conditions for both surgeon and OR-assistant, stable image and the reduction of the number of people needed for a MIS-procedure.

## Product:

Mofix fixates endoscopes (and other instruments) in any desired position. By pushing a switch on the endoscope it can be repositioned single handed and in an intuitive way, "drag and drop". The patented Mofix technology is currently successfully in use at various types of laparoscopic surgery, i.e. urology, gynecology and general surgery.



[mofix.com](http://mofix.com)

### Mofix bv

Pantheon 26  
7521PR Enschede  
Netherlands

T +31(0)534803927  
M +31(0)651577276  
F +31(0)534803922  
info@mofix.com

### Distributor Netherlands

Van Straten Medical BV  
Edisonbaan 20  
3430AK Nieuwegein, NL

T +31(0)30 6023830  
vanstratenmedical.com  
oktechniek@vanstratenmedical.com